

BEFORE / AFTER | 5 months post 7 tx | courtesy of Carrie Labert, PA-C



Forever Clear BBL™

Step 1: 420 nm Filter, 5 J/cm², 260 ms, 15 °C Cooling, 3 Passes

Step 2: 560 nm Filter, 9-12 J/cm², 15-20 ms, 15 °C Cooling, 2 Passes

560 nm Filter, 15 J/cm², 18 ms, 15 °C Cooling, 2-3 Passes (Checks, Nose, Chin)

Step 3: SkinTye: 590 Filter, 12 W/cm², 12 sec, 30 °C Cooling

