

BEFORE / AFTER | 4 months post 6 tx | courtesy of Patrick Bitter, Jr., MD



# Forever Clear BBL™

Step 1: 420 nm filter, 5 J/cm<sup>2</sup>, 240 ms, 15 °C cooling, 3 passes with large spot

Step 2: 560 nm filter, 15 J/cm<sup>2</sup>, 30 ms, 15 °C cooling, 2 Passes with square adapter

Step 3: 590 ST filter, 10 W/cm<sup>2</sup>, 12 s, 30 °C cooling, 8 passes with large spot

